

Dilate salivary duct strictures.




Illustration by Lisa Clark

Advance[®]

SALIVARY DUCT BALLOON CATHETER SET

High-strength

The high-strength, minimally compliant material allows for the complete dilation of the stricture.

Minimally invasive

The flexible, low-profile design supports minimally invasive intervention through radial dilation.

Multi-modal visualization

The catheter is compatible with direct endoscopic visualization or ultrasound imaging based on physician preference.





Advance®

SALIVARY DUCT BALLOON CATHETER SET

Used for the dilation of strictures in the submandibular and parotid salivary ducts.

The balloon is constructed of an extra-thin wall and high-strength, minimally compliant material that is intended to be used with sialendoscopic visualization in the submandibular salivary ducts or with ultrasound imaging in parotid or submandibular salivary ducts.

Components

- Advance Salivary Duct Balloon Catheter
- Wire guide

Order Number	Reference Part Number	Wire Guide Diameter inch	Inflated Balloon Diameter mm	Balloon Length cm	Preinflated Balloon Diameter mm
G34991	SDBS-14-30-1.5-2	.014	1.5	2	1
G34992	SDBS-14-30-2.0-2	.014	2.0	2	1
G34993	SDBS-14-30-2.5-2	.014	2.5	2	1
G34994	SDBS-14-30-3.0-2	.014	3.0	2	1

Related products



Cook Sphere® Inflation Device

Order Number	Reference Part Number	Fluid Capacity mL	Maximum Pressure atm
G34903	CID-60-15	60	15

Some products or part numbers may not be available in all markets. Contact your local Cook Medical sales representative or Customer Service for details.

Customer Service

EU Website: cookmedical.eu
 EDI: cookmedical.eu/edi
 Distributors: +353 61239240, ssc.distributors@cookmedical.com
 Austria: +43 179567121, oe.orders@cookmedical.com
 Belgium: +32 27001702, be.orders@cookmedical.com
 Denmark: +45 38487607, da.orders@cookmedical.com
 Finland: +358 972519996, fi.orders@cookmedical.com
 France: +33 171230269, fr.orders@cookmedical.com
 Germany: +49 6950072804, de.orders@cookmedical.com
 Hungary: +36 17779199, hu.orders@cookmedical.com
 Iceland: +354 800 7615, IS.orders@cookmedical.com
 Ireland: +353 61239252, ie.orders@cookmedical.com
 Italy: +39 0269682853, it.orders@cookmedical.com
 Netherlands: +31 202013367, nl.orders@cookmedical.com
 Norway: +47 23162968, no.orders@cookmedical.com
 Spain: +34 912702691, es.orders@cookmedical.com
 Sweden: +46 858769468, se.orders@cookmedical.com
 Switzerland - French: +41 448009609, fr.orders@cookmedical.com
 Switzerland - Italian: +41 448009609, it.orders@cookmedical.com
 Switzerland - German: +41 448009609, de.orders@cookmedical.com
 United Kingdom: +44 2073654183, uk.orders@cookmedical.com

USA Website: cookmedical.com
 EDI: cookmedical.com/edi.do
Americas:
 Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335
 E-mail: orders@cookmedical.com
Australia:
 Phone: +61 734346000, 1800777222, Fax: +61 734346001, 1800077283
 E-mail: cau.custserv@cookmedical.com



Exclusive UK distributor:

Est: 1966 - supplying NHS & MoD

Tel: 01264 332172 • E-Mail: info@ccmed.co.uk
 Web: www.CCMed.co.uk • Twitter: [@CCmedUK](https://twitter.com/CCmedUK)

AI, ESC, IR, OHNS, PI, RH, SUR-A4